



 **Ambix® nova**

Designed for life

**User guide**

Ambulatory pump, giving set and accessories

 **Calea**  
the heart of homecare

 **FRESENIUS  
KABI**  
caring for life

# Welcome to Ambix nova

Ambix nova has been chosen by your team of healthcare professionals as a user-friendly solution that provides optimal safety and simple infusion at home.

While your Homecare nurse will have explained the relevant information you need for your specific treatment and pump, this booklet supports those instructions and your use of the Ambix nova pump. Please refer to the instructions for use (IFU) for full instructions.

## What will be in the box

- Pump
- Smart holder power
- Mains lead
- Instructions for use
- Quick reference guide

## Before connecting



Charge the battery for at least six hours: pump has a minimum battery life of 24 hours at 25 mL/h














Prime the line, as you were shown in your nurse training

## Servicing and support

If you have any issues whilst using your Ambix nova pump, you can use the troubleshooting guide within this booklet. If you cannot solve your issue, please call the Freephone 24-hour advice line on **0800 121 8300**.


Please see the back of your pump for the date of your next service. The Calea Patient Services team will be in touch a few weeks before this date to arrange the pump service.

# Ambix nova pump


- Alarm silence 
- Keypad lock 
- Menu 
- Mains supply (LED icon below word 'nova' which goes green when connected)
- Switch the pump on 
- Jump to another parameter/next option  
- Adjust parameter values   keys
- Start the pump 
- Stop the pump or silence the alarm 
- Priming 




## Pump Back




## Pump Side



## Smart holder power Front



## Smart holder power Back



# Ambix nova giving set





# How to install & remove your Ambix nova giving set

## Installing the set









- Open the door with the lever
- Match the green part of the Ambix nova giving set to the green part of the pump - make sure the foot of the clamp is facing outwards
- Match the blue top of the anti-kink tubing on the Ambix nova giving set to the blue part of the pump
- Close the door and push the lever back into place

## Removing the set

- Stop the pump by pressing STOP
- To turn off the pump you need to press ON/OFF for at least 3 seconds
- The pump will count down 3...2...1...off
- Open the pump door
- Remove the Ambix nova giving set
- Dispose of the set following the correct procedure, as advised in your nurse training





# How to use your Ambix nova pump

- Switch the pump on using the ON/OFF button 
- Press MENU  to choose the correct mode (please see relevant mode pages for more information)
- Check that the settings match your prescription or change to make sure they match
- Start the infusion using the START button 
- Stop the pump using the STOP button 
- Press the MENU  button to retrieve the last volume of fluids infused
- Press the MENU  button and hold for 3 seconds to zero the volume infused, if required
- Use the alarm silence button or press STOP  to silence the alarm 



To select your infusion mode, press MENU  before commencing infusion



Select the infusion mode icon, which looks like this:  and validate by pressing the green button  to enter the infusion mode selection menu











Select the required   infusion mode and press the green button 

# How to set your specific mode






## Flowrate only

**Allows only one setting to be adjusted - flowrate**

1. Select flowrate only mode: .  
Use the arrow keys   to select and green button  to OK
2. The flowrate can be adjusted using the   keys
3. Press START  to start the infusion
4. Monitoring screen will show infused volume
5. Press MENU  for more infusion information








## Continuous

**Enables continuous infusion at a pre-set flow rate, volume and duration**

1. Select flowrate only mode: .  
Use the arrow keys to select and green button to OK
2. The flowrate can be adjusted using the   keys
3. Press START  to start the infusion
4. Monitoring screen will show infused volume
5. Press MENU  for more infusion information

## Ramp

**Enables continuous infusion at a pre-set flow rate, volume and duration**

1. Select ramp mode: .  
Use the arrow keys to select and green button to OK
2. Programming screen will show volume, ramp up duration, ramp down duration and total duration
3. Adjust according to prescription volume, ramp up duration, ramp down duration and total duration using   keys and arrow keys   to move between settings. Flow rate will be calculated by the pump
4. Press START  to start the infusion
5. Monitoring screen will show phase, remaining time, infused volume and total volume to be infused
6. Press MENU  for more infusion information



# Alarms & safety features

The Ambix nova pump will provide alarm and pre-alarm information, when needed.

When an alarm occurs, the pump will show a symbol and message and the top left-hand light will flash to indicate the reason for the alarm.

An alarm will also sound - you can press the alarm symbol button to mute the beeping or press STOP to release the alarm and clear the signal.

## Lights



**Green light** for pumping indication (not alarm)




**Yellow light** for low priority alarm



**Red light** for high priority alarm

## Alarm system


### Low priority



#### Empty battery pre-alarm:

Battery charge has 30 minutes remaining.

**Remedy:** Plug the pump into the mains.



#### End of infusion pre-alarm:

Target volume will soon be reached.

**Remedy:** Stop the infusion when finished and clean the pump.




## High priority




**Missing set alarm:** Missing set, set not properly installed or wrong set installed.

**Remedy:** Check and/or reinstall set. If alarm continues, contact your Homecare company.



**Door open alarm:** Pump door not properly closed or door is removed or faulty.

**Remedy:** Close the pump door and the lever. Remove the door and re-hang. Contact your Homecare company.




**Up-stream occlusion:** Upstream flow path is blocked between the container and the pump.

**Remedy:** Remove the kink, blockages, clamps or air particles.




**Down-stream occlusion:** Downstream flow path is blocked below the pump at the patient side.

**Remedy:** Remove the kink, blockages, clamps or air particles.




**Empty bag / air in line alarm:** Infusion container is empty or air in the set.

**Remedy:** Stop the infusion / change / check the bag and Ambix giving set. Prime the Ambix giving set. Clean the pump.



**End of infusion alarm:** Infusion is finished and the target volume is reached.

**Remedy:** After the infusion has finished, clean the pump.



**Empty battery alarm:** Battery is almost empty.

**Remedy:** Plug the pump into the mains.



## Technical alarm

**Remedy:** Switch the pump off and contact your Homecare company.



## Accessing history menu

1. Press MENU
2. Choose History: Use the arrow keys to select and green button to OK
3. Press Start/OK button
4. Choose the item you need using the arrow keys to select and green button to OK
5. Press Start/OK button to view the information

# How to clean your pump

It is important to clean your pump:



Before connecting  
it to the mains



When fluids or drugs come  
into contact with the pump



At least once a week using  
the recommended cleaning  
and disinfection products  
(detailed in the IFU)

## Instructions for cleaning

- **STEP 1**  
Thoroughly wipe down all surfaces of the pump from top to bottom
- **STEP 2**  
Gently wipe down the mechanism and sensor area using a ready-to-use wipe
- **STEP 3**  
Use cotton wool to clean contact pins
- **STEP 4**  
Remove the pump door then repeat Step 1 with the pump door and holder
- **STEP 5**  
Wipe down the power cord
- **STEP 6**  
Allow the device to dry completely at room temperature. Re-install the pump door before use



## WARNING

- Do not submerge the pump in water or any other liquid. If this happens, contact your homecare company
- Do not use in the dishwasher
- Do not sterilise the pump



# Ambix nova: a complete solution

**Battery life:**  
minimum of 24 hours at  
25 mL/h

**Charging:**  
6 hours with pump off

**Flow rates:**  
Min 1 mL/h and  
Max 600 mL/h

**Power consumption:**  
maximum 9W in  
standard operating  
conditions

**Safe to take through  
security scanners**



Calea, Cestrian Court, Eastgate Way, Manor Park, Runcorn, Cheshire, WA7 1NT  
Tel: 01928 533500 Freephone 24-hour advice line 0800 121 8300

**[www.calea.co.uk](http://www.calea.co.uk)**

© Fresenius Kabi Limited. UK-CAL-2100070 | September 2022